



Endoscope Washer OLYMPUS MINI ETD (Reconditioned)

- Condition: used - reconditioned,
- Made in Germany,
- Automatic,
- 400 V 3 phase power supply,
- Frequency 50-60 Hz,
- Chamber for 1 endoscope,
- Built-in digital display,
- Solid construction made of stainless steel,
- Work in a closed system - a completely sealed structure excludes harmful fumes originating from disinfectants,
- Connections for Olympus endoscopes,
- Automatic endoscope disinfection process,
- Chamber dimensions: 39 x 62 x 48 cm,
- Dimensions: 60 x 63 x 84 cm,
- Technical Passport,
- 30-days warranty

