



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Anamed ECG Cable (New)

- New,
- Made in Uruguay,
- Code: IEC
- Options:
 - **Cable Length: 3,5 m**
 - 3 leads snap (BIOECG.BC.S3F),
 - 3 leads clamp (BIOECG.BC.C3F),
 - 5 leads snap (BIOECG.BC.S5F),
 - 5 leads clamp (BIOECG.BC.C5F),
 - **Compatible with:**
 - AM-78100A, AM 78100B, AM 78600, Vital Lite, Vital Wave, AM-7810A,
 - **Cable length: 4,0 m**
 - 3 leads snap (BIOECG.BC.S3D),
 - 3 leads clamp (BIOECG.BC.C3D),
 - 5 leads snap (BIOECG.BC.S5D),
 - 5 leads clamp (BIOECG.BC.C5D),
 - **Compatible with:**
 - AM-78100A, AM 78100B, AM 78600, Vital Lite, Vital Wave, AM-7810A,
 - **Cable Length: 3,0 m (this product requires a lead set)**
 - 3 leads (BIOECG.BC.3T),
 - 5 leads (BIOECG.BC.5T),
 - **Compatible with:**
 - AM-78100A, AM 78100B, AM 78600, Vital Lite, Vital Wave, AM-7810A,
 - **Cable Length: 1,0 m (this product requires a trunk cable)**
 - 3 leads snap (BIOECG.RA.S3),
 - 3 leads clamp (BIOECG.RA.C3),
 - 5 leads snap (BIOECG.RA.S5),
 - 5 leads clamp (BIOECG.RA.C5),
 - **Compatible with:**
 - AM-78100A, AM 78100B, AM 78600, Vital Lite, Vital Wave, AM-7810A,
- **Cable extension available on demand,**
- 12 months warranty



Accessories	Options	Price
ECG lead set 3 or 5 leads	Lead set	2000000000
ECG lead set 3 or 5 leads	Lead set	2000000000
ECG lead set 3 or 5 leads	Lead set	2000000000
ECG lead set 3 or 5 leads	Lead set	2000000000

Accessories	Options	Price
ECG lead set 3 or 5 leads	Lead set	2000000000
ECG lead set 3 or 5 leads	Lead set	2000000000
ECG lead set 3 or 5 leads	Lead set	2000000000
ECG lead set 3 or 5 leads	Lead set	2000000000

Accessories	Options	Price
ECG lead set 3 or 5 leads	Lead set	2000000000
ECG lead set 3 or 5 leads	Lead set	2000000000
ECG lead set 3 or 5 leads	Lead set	2000000000
ECG lead set 3 or 5 leads	Lead set	2000000000

Accessories	Options	Price
ECG lead set 3 or 5 leads	Lead set	2000000000
ECG lead set 3 or 5 leads	Lead set	2000000000
ECG lead set 3 or 5 leads	Lead set	2000000000
ECG lead set 3 or 5 leads	Lead set	2000000000