



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Burdick ECG Cable (New)

- New,
- Made in Uruguay,
- Code: IEC
- Options:
 - **Cable Length: 3,0 m**
 - 10 lead Banana 4 mm (BIOECG.BU.P10F),
 - **Compatible with:**
 - EK10, Elite, E350, E350I, E550, E560, Eclipse LE, LE I, Elite II, Quinton Eclipse Plus,
 - **Cable Length: 3,0 m**
 - 10 lead Banana 4 mm (BIOECG.BU2.P10F),
 - **Compatible with:**
 - Eclipse 4, 300, 400, 800, 850, Plus and Premier, Elite II, Eclipse LE II and Atria Systems, Atria 3000, 3100, 6100,
- **Cable extension available on demand,**
- 12 months warranty

