



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Cardioline ECG Cable (New)

- New,
- Made in Uruguay,
- Code: IEC
- Options:
 - **Cable Length: 0,6 m**
 - 5 leads snap (BIOECG.CDK.S5F),
 - 5 leads clamp (BIOECG.SDK.S5F),
 - **Compatible with:**
 - Click Holter,
 - **Cable Length: 0,6 m**
 - 10 lead banana (BIOECG.CAK.P10F)
 - **Compatible with:**
 - Click Holter,
 - **Cable Length: 3,0 m**
 - 10 lead banana (BIOECG.NK.P10F)
 - **Compatible with:**
 - E1, ER1, E3, ITA300, C3,
 - **Cable Length: 3,0 m**
 - 10 lead banana (BIOECG.CA.P10F),
 - **Compatible with:**
 - ITA 1 -TA1 Plus, ITA 3 - ITA 3 Plus, ITA 30, ITA 60 - ITA 60 Plus, Elan, Elan 100, Elan 1100, 103, 106, 107, 302, 303, 1201 and Cardioline Remco Delta 1, Delta 3, Delta 30, Delta 30D, Delta 1 Plus, Delta 3 Plus, Delta 60 Plus, ECG 106,
- **Cable extension available on demand,**
- 12 months warranty

