



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Datascope ECG Cable (New)

- New,
- Made in Uruguay,
- Code: IEC
- Options:
 - **Cable Length: 3,5 m**
 - 3 leads snap (BIOECG.BC.S3F),
 - 3 leads clamp (BIOECG.BC.C3F),
 - 5 leads snap (BIOECG.BC.S5F),
 - 5 leads clamp (BIOECG.BC.C5F),
 - **Compatible with:**
 - Datascope 2000, 2000A, 2001, 2002, 2100, 2200, 3000, 3000A, Ballon Pump, M/D3, M/D3A, M/D4, Passport, Passport IL, Passport EL, Passport XG, Passport XL, Passport 2, Passport 2LT, Passport 5L, Trio, Spectrum, Spectrum OR, System 82, System 83, System 90, System 90T, System 96, System 97,
 - **Cable length: 4,0 m**
 - 3 leads snap (BIOECG.BC.S3D),
 - 3 leads clamp (BIOECG.BC.C3D),
 - 5 leads snap (BIOECG.BC.S5D),
 - 5 leads clamp (BIOECG.BC.C5D),
 - **Compatible with:**
 - Datascope 2000, 2000A, 2001, 2002, 2100, 2200, 3000, 3000A, Ballon Pump, M/D3, M/D3A, M/D4, Passport, Passport IL, Passport EL, Passport XG, Passport XL, Passport 2, Passport 2LT, Passport 5L, Trio, Spectrum, Spectrum OR, System 82, System 83, System 90, System 90T, System 96, System 97,
 - **Cable Length: 3,0 m (this product requires a lead set)**
 - 3 leads (BIOECG.BC.3T),
 - 5 leads (BIOECG.BC.5T),
 - **Compatible with:**
 - Datascope 2000, 2000A, 2001, 2002, 2100, 2200, 3000, 3000A, Ballon Pump, M/D3, M/D3A, M/D4, Passport, Passport IL, Passport EL, Passport XG, Passport XL, Passport 2, Passport 2LT, Passport 5L, Trio, Spectrum, Spectrum OR, System 82, System 83, System 90, System 90T, System 96, System 97,
 - **Cable Length: 1,0 m (this product requires a trunk cable)**
 - 3 leads snap (BIOECG.RA.S3),
 - 3 leads clamp (BIOECG.RA.C3),
 - 5 leads snap (BIOECG.RA.S5),
 - 5 leads clamp (BIOECG.RA.C5),
 - **Compatible with:**
 - Datascope 2000, 2000A, 2001, 2002, 2100, 2200, 3000, 3000A, Ballon Pump, M/D3, M/D3A, M/D4, Passport, Passport IL, Passport EL, Passport XG, Passport XL, Passport 2, Passport 2LT, Passport 5L, Trio, Spectrum, Spectrum OR, System 82, System 83, System 90, System 90T, System 96, System 97,
- **Cable extension available on demand,**
- 12 months warranty

