



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



EKG-Edan cable (New)

New,
Uruguay production,
Number of leads: 10,
Encoding color: IEC
Plug type: 4 mm banana,
Cable length 3.0 m (with the option of ordering a cable of any length),
The cable is used in series devices:
All devices from Edan,
12-month warranty

