



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordynator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



### Elettronica Trentina-Cardiette ECG Cable (New)

- New,
- Made in Uruguay,
- 10 lead cable,
- Code: IEC
- Plug: banana 4 mm,
- **Cable Length 3,0 m** (cable extension available on demand),
- Compatible with:
  - All Edan devices,
- 12 months warranty

