



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Email-Transmai ECG Cable (New)

- New,
- Made in Uruguay,
- Code: IEC
- Options:
 - **Cable Length: 3,5 m**
 - 3 leads snap (BIOECG.BC.S3F),
 - 3 leads clamp (BIOECG.BC.C3F),
 - 5 leads snap (BIOECG.BC.S5F),
 - 5 leads clamp (BIOECG.BC.C5F),
 - **Compatible with:**
 - RX 300 y MX series,
 - **Cable length: 4,0 m**
 - 3 leads snap (BIOECG.BC.S3D),
 - 3 leads clamp (BIOECG.BC.C3D),
 - 5 leads snap (BIOECG.BC.S5D),
 - 5 leads clamp (BIOECG.BC.C5D),
 - **Compatible with:**
 - RX 300 y MX series,
 - **Cable Length: 3,0 m (this product requires a lead set)**
 - 3 leads (BIOECG.BC.3T),
 - 5 leads (BIOECG.BC.5T),
 - **Compatible with:**
 - RX 300 y MX series,
 - **Cable Length: 1,0 m (this product requires a trunk cable)**
 - 3 leads snap (BIOECG.RA.S3),
 - 3 leads clamp (BIOECG.RA.C3),
 - 5 leads snap (BIOECG.RA.S5),
 - 5 leads clamp (BIOECG.RA.C5),
 - **Compatible with:**
 - RX 300 y MX series,
- Cable extension available on demand,
- 12 months warranty



Email / Transmai	Options	Code
3 leads snap	3,5 m	BIOECG.BC.S3F
3 leads clamp	3,5 m	BIOECG.BC.C3F
5 leads snap	3,5 m	BIOECG.BC.S5F
5 leads clamp	3,5 m	BIOECG.BC.C5F

Email / Transmai	Options	Code
3 leads snap	4,0 m	BIOECG.BC.S3D
3 leads clamp	4,0 m	BIOECG.BC.C3D
5 leads snap	4,0 m	BIOECG.BC.S5D
5 leads clamp	4,0 m	BIOECG.BC.C5D

Email / Transmai	Options	Code
3 leads snap	3,0 m	BIOECG.BC.3T
3 leads clamp	3,0 m	BIOECG.BC.C3T
5 leads snap	3,0 m	BIOECG.BC.5T
5 leads clamp	3,0 m	BIOECG.BC.C5T

Email / Transmai	Options	Code
3 leads snap	1,0 m	BIOECG.RA.S3
3 leads clamp	1,0 m	BIOECG.RA.C3
5 leads snap	1,0 m	BIOECG.RA.S5
5 leads clamp	1,0 m	BIOECG.RA.C5