



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Esaote ECG Cable (New)

- New,
- Made in Uruguay,
- Code: IEC
- Options:
 - **Cable Length: 3,5 m**
 - 3 leads snap (BIOECG.ES.S3F),
 - 3 leads clamp (BIOECG.ES.C3F),
 - **Compatible with:**
 - All Esaote devices,
 - **Cable Length: 3,0 m**
 - 10 lead banana (BIOECG.SC.P10F),
 - **Compatible with:**
 - EKG P120,
 - **Cable Length: 3,0 m**
 - 10 lead banana (BIOECG.SC2.P10F),
 - **Compatible with:**
 - EKG P80, P8000 Series,
- **Cable extension available on demand,**
- 12 months warranty

