



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Feas ECG Cable (New)

- New,
- Made in Uruguay,
- Code: IEC
- Options:
 - **Cable Length: 3,5 m**
 - 3 leads snap (BIOECG.FE.S3F),
 - 3 leads clamp (BIOECG.FE.C3F),
 - 5 leads snap (BIOECG.FE.S5F),
 - 5 leads clamp (BIOECG.FE.C5F),
 - **Compatible with:**
 - All Feas devices,
 - **Cable length: 4,0 m**
 - 3 leads snap (BIOECG.FE.S3D),
 - 3 leads clamp (BIOECG.FE.C3D),
 - 5 leads snap (BIOECG.FE.S5D),
 - 5 leads clamp (BIOECG.FE.C5D),
 - **Compatible with:**
 - All Feas devices,
 - **Cable Length: 3,0 m (this product requires a lead set)**
 - 3 leads (BIOECG.FE.3T),
 - 5 leads (BIOECG.FE.5T),
 - **Compatible with:**
 - All Feas devices,
 - **Cable Length: 1,0 m (this product requires a trunk cable)**
 - 3 leads snap (BIOECG.RA.S3),
 - 3 leads clamp (BIOECG.RA.C3),
 - 5 leads snap (BIOECG.RA.S5),
 - 5 leads clamp (BIOECG.RA.C5),
 - **Compatible with:**
 - All Feas devices,
- Cable extension available on demand,
- 12 months warranty



Feas	Options	Code
3/5 leads snap with 3/5 m 1.0 m cable	3 leads snap	BIOECG.FE.S3F
	5 leads snap	BIOECG.FE.S5F
3/5 leads clamp with 3/5 m 1.0 m cable	3 leads clamp	BIOECG.FE.C3F
	5 leads clamp	BIOECG.FE.C5F

Feas	Options	Code
3/5 leads snap with 3/5 m 1.0 m cable	3 leads snap	BIOECG.FE.S3D
	5 leads snap	BIOECG.FE.S5D
3/5 leads clamp with 3/5 m 1.0 m cable	3 leads clamp	BIOECG.FE.C3D
	5 leads clamp	BIOECG.FE.C5D

Feas	Options	Code
3/5 leads snap with 3/5 m 1.0 m cable	3 leads snap	BIOECG.FE.3T
	5 leads snap	BIOECG.FE.5T

Feas	Options	Code
3/5 leads snap with 3/5 m 1.0 m cable	3 leads snap	BIOECG.RA.S3
	5 leads snap	BIOECG.RA.S5
3/5 leads clamp with 3/5 m 1.0 m cable	3 leads clamp	BIOECG.RA.C3
	5 leads clamp	BIOECG.RA.C5