



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordynator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **Fukuda Denshi ECG Cable (New)**

- New,
- Made in Uruguay,
- 10 lead cable,
- Code: IEC
- Plug: banana 4 mm,
- **Cable Length 3,0 m** (cable extension available on demand),
- Compatible with:
  - FX-101, FX-102F, FX-3010, FX-328UY, FCP 2155, FCP 2201, FCP 4101, FCP 4101A, CardiMax FX-7101,FX-7202, FX-7402,
- 12 months warranty

