



PRAISTON Sp. z o.o.
ul. Górska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybierański
+48 882 762 006

kw@praiston.pl



GE-Datex Ohmeda ECG Cable (New)

- New,
- Made in Uruguay,
- Code: IEC
- Options:
 - **Cable Length: 3,5 m**
 - 3 leads snap (BIOECG.DA.S3F),
 - 3 leads clamp (BIOECG.DA.C3F),
 - 5 leads snap (BIOECG.DA.S5F),
 - 5 leads clamp (BIOECG.DA.C5F),
 - **Compatible with:**
 - AS/3, CS/3, S/5, Cardiocap I, Cardiocap II, Cardiocap 5, Light Monitor, Satlite plus,
 - **Cable Length: 4,0 m**
 - 3 leads snap (BIOECG.DA.S3D),
 - 3 leads clamp (BIOECG.DA.C3D),
 - 5 leads snap (BIOECG.DA.S5D),
 - 5 leads clamp (BIOECG.DA.C5D),
 - **Compatible with:**
 - AS/3, CS/3, S/5, Cardiocap I, Cardiocap II, Cardiocap 5, Light Monitor, Satlite plus,
 - **Cable Length: 3,0 m (this product requires a lead set)**
 - 3 leads (BIOECG.DA.3T),
 - 5 leads (BIOECG.DA.5T),
 - **Compatible with:**
 - AS/3, CS/3, S/5, Cardiocap I, Cardiocap II, Cardiocap 5, Light Monitor, Satlite plus,
 - **Cable Length: 1,0 m (this product requires a trunk cable)**
 - 3 leads snap (BIOECG.RA.S3),
 - 3 leads clamp (BIOECG.RA.C3),
 - 5 leads snap (BIOECG.RA.S5),
 - 5 leads clamp (BIOECG.RA.C5),
 - **Compatible with:**
 - AS/3, CS/3, S/5, Cardiocap I, Cardiocap II, Cardiocap 5, Light Monitor, Satlite plus,
- **Cable extension available on demand,**
- 12 months warranty

