



PRAISTON Sp. z o.o.
ul. Górska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



GE-Marquette ECG Cable (New)

- New,
- Made in Uruguay,
- Code: IEC
- Options:
 - **Cable Length: 3,5 m**
 - 3 leads snap (BIOECG.MQ.S3F),
 - 3 leads clamp (BIOECG.MQ.C3F),
 - 5 leads snap (BIOECG.MQ.S5F),
 - 5 leads clamp (BIOECG.MQ.C5F),
 - **Compatible with:**
 - Dash 1000, Dash Pro 2000, Dash 2500, Dash Pro 3000, Dash Pro 4000, Dash 5000, Eagle Systems, MACLab, Solar 3000, Solar 8000, Solar 9500, Aestiva/5, Tram 100, Tram 200, Tram 300, Tram 400, Tram 500, Tram 600, Tram 800, Unity Systems and Carescape,
 - **Cable length: 4,0 m**
 - 3 leads snap (BIOECG.MQ.S3D),
 - 3 leads clamp (BIOECG.MQ.C3D),
 - 5 leads snap (BIOECG.MQ.S5D),
 - 5 leads clamp (BIOECG.MQ.C5D),
 - **Compatible with:**
 - Dash 1000, Dash Pro 2000, Dash 2500, Dash Pro 3000, Dash Pro 4000, Dash 5000, Eagle Systems, MACLab, Solar 3000, Solar 8000, Solar 9500, Aestiva/5, Tram 100, Tram 200, Tram 300, Tram 400, Tram 500, Tram 600, Tram 800, Unity Systems and Carescape,
 - **Cable Length: 3,0 m (this product requires a lead set)**
 - 3 leads (BIOECG.MQ.3T),
 - 5 leads (BIOECG.MQ.5T),
 - **Compatible with:**
 - Dash 1000, Dash Pro 2000, Dash 2500, Dash Pro 3000, Dash Pro 4000, Dash 5000, Eagle Systems, MACLab, Solar 3000, Solar 8000, Solar 9500, Aestiva/5, Tram 100, Tram 200, Tram 300, Tram 400, Tram 500, Tram 600, Tram 800, Unity Systems and Carescape,
 - **Cable Length: 1,0 m (this product requires a trunk cable)**
 - 3 leads snap (BIOECG.RA.S3),
 - 3 leads clamp (BIOECG.RA.C3),
 - 5 leads snap (BIOECG.RA.S5),
 - 5 leads clamp (BIOECG.RA.C5),
 - **Compatible with:**
 - Dash 1000, Dash Pro 2000, Dash 2500, Dash Pro 3000, Dash Pro 4000, Dash 5000, Eagle Systems, MACLab, Solar 3000, Solar 8000, Solar 9500, Aestiva/5, Tram 100, Tram 200, Tram 300, Tram 400, Tram 500, Tram 600, Tram 800, Unity Systems and Carescape,
 - **Cable Length: 3,0 m**
 - 10 lead banana 4 mm (BIOECG.MQ.P10F),
 - **Compatible with:**
 - MAC 400, MAC 500, MAC 1100, MAC 1200,
- **Cable extension available on demand,**
- 12 months warranty

