



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Mennen ECG Cable (New)

- New,
- Made in Uruguay,
- Code: IEC
- Options:
 - **Cable Length: 3,5 m**
 - 3 leads snap (BIOECG.ME.S3F),
 - 3 leads clamp (BIOECG.ME.C3F),
 - 5 leads snap (BIOECG.ME.S5F),
 - 5 leads clamp (BIOECG.ME.C5F),
 - **Compatible with:**
 - Horizon 250, 450, 740, 741, 742, 3000, Microsolo, 936 Series, 425/101 Preamp, Cardio-Pak 2000,
 - **Cable Length: 4,0 m**
 - 3 leads snap (BIOECG.ME.S3D),
 - 3 leads clamp (BIOECG.ME.C3D),
 - 5 leads snap (BIOECG.ME.S5D),
 - 5 leads clamp (BIOECG.ME.C5D),
 - **Compatible with**
 - Horizon 250, 450, 740, 741, 742, 3000, Microsolo, 936 Series, 425/101 Preamp, Cardio-Pak 2000,
 - **Cable Length: 3,0 m (this product requires a lead set)**
 - 3 leads (BIOECG.ME.3T),
 - 5 leads (BIOECG.ME.5T),
 - **Compatible with:**
 - Horizon 250, 450, 740, 741, 742, 3000, Microsolo, 936 Series, 425/101 Preamp, Cardio-Pak 2000,
 - **Cable Length: 1,0 m (this product requires a trunk cable)**
 - 3 leads snap (BIOECG.RA.S3),
 - 3 leads clamp (BIOECG.RA.C3),
 - 5 leads snap (BIOECG.RA.S5),
 - 5 leads clamp (BIOECG.RA.C5),
 - **Compatible with:**
 - Horizon 250, 450, 740, 741, 742, 3000, Microsolo, 936 Series, 425/101 Preamp, Cardio-Pak 2000,
 - **Cable Length: 3,5 m**
 - 3 leads snap (BIOECG.ME3.S3F),
 - 3 leads clamp (BIOECG.ME3.C3F),
 - 5 leads snap (BIOECG.ME3.S5F),
 - 5 leads clamp (BIOECG.ME3.C5F),
 - **Compatible with:**
 - Enmove, 1000, 1200, Envoy, Vitalogik,
 - **Cable Length: 4,0 m**
 - 3 leads snap (BIOECG.ME3.S3D),
 - 3 leads clamp (BIOECG.ME3.C3D),
 - 5 leads snap (BIOECG.ME3.S5D),
 - 5 leads clamp (BIOECG.ME3.C5D),
 - **Compatible with:**
 - Enmove, 1000, 1200, Envoy, Vitalogik,
 - **Cable Length: 3,0 m (this product requires a lead set)**
 - 3 leads (BIOECG.ME3.3T),
 - 5 leads (BIOECG.ME3.5T),
 - **Compatible with:**
 - Enmove, 1000, 1200, Envoy, Vitalogik,
 - **Cable Length: 1,0 m (this product requires a trunk cable)**
 - 3 leads snap (BIOECG.RA.S3),
 - 3 leads clamp (BIOECG.RA.C3),
 - 5 leads snap (BIOECG.RA.S5),
 - 5 leads clamp (BIOECG.RA.C5),
 - **Compatible with:**
 - Enmove, 1000, 1200, Envoy, Vitalogik,
- **Cable extension available on demand,**
- 12 months warranty



Hansen		
Options	Price	Code
BNC connector cable	1.00m	8000040000
2 x 1.00m	1.00m	8000040000
2 x 1.00m	1.00m	8000040000
2 x 1.00m	1.00m	8000040000
2 x 1.00m	1.00m	8000040000

Hansen		
Options	Price	Code
BNC connector cable	1.00m	8000040000
2 x 1.00m	1.00m	8000040000
2 x 1.00m	1.00m	8000040000
2 x 1.00m	1.00m	8000040000
2 x 1.00m	1.00m	8000040000

Hansen		
Options	Price	Code
BNC connector cable	1.00m	8000040000
2 x 1.00m	1.00m	8000040000
2 x 1.00m	1.00m	8000040000
2 x 1.00m	1.00m	8000040000
2 x 1.00m	1.00m	8000040000

Hansen		
Options	Price	Code
BNC connector cable	1.00m	8000040000
2 x 1.00m	1.00m	8000040000
2 x 1.00m	1.00m	8000040000
2 x 1.00m	1.00m	8000040000
2 x 1.00m	1.00m	8000040000

Hansen		
Options	Price	Code
BNC connector cable	1.00m	8000040000
2 x 1.00m	1.00m	8000040000
2 x 1.00m	1.00m	8000040000
2 x 1.00m	1.00m	8000040000
2 x 1.00m	1.00m	8000040000

Hansen		
Options	Price	Code
BNC connector cable	1.00m	8000040000
2 x 1.00m	1.00m	8000040000
2 x 1.00m	1.00m	8000040000
2 x 1.00m	1.00m	8000040000
2 x 1.00m	1.00m	8000040000

Hansen		
Options	Price	Code
BNC connector cable	1.00m	8000040000
2 x 1.00m	1.00m	8000040000
2 x 1.00m	1.00m	8000040000
2 x 1.00m	1.00m	8000040000
2 x 1.00m	1.00m	8000040000

Hansen		
Options	Price	Code
BNC connector cable	1.00m	8000040000
2 x 1.00m	1.00m	8000040000
2 x 1.00m	1.00m	8000040000
2 x 1.00m	1.00m	8000040000
2 x 1.00m	1.00m	8000040000