



PRAISTON Sp. z o.o.
ul. Górska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Moriya ECG Cable (New)

- New,
- Made in Uruguay,
- Code: IEC
- Options:
 - **Cable Length: 3,5 m**
 - 3 leads snap (BIOECG.BC.S3F),
 - 3 leads clamp (BIOECG.BC.C3F),
 - 5 leads snap (BIOECG.BC.S5F),
 - 5 leads clamp (BIOECG.BC.C5F),
 - **Compatible with:**
 - M2000, M3000, M9100,
 - **Cable length: 4,0 m**
 - 3 leads snap (BIOECG.BC.S3D),
 - 3 leads clamp (BIOECG.BC.C3D),
 - 5 leads snap (BIOECG.BC.S5D),
 - 5 leads clamp (BIOECG.BC.C5D),
 - **Compatible with:**
 - M2000, M3000, M9100,
 - **Cable Length: 3,0 m (this product requires a lead set)**
 - 3 leads (BIOECG.BC.3T),
 - 5 leads (BIOECG.BC.5T),
 - **Compatible with:**
 - M2000, M3000, M9100,
 - **Cable Length: 1,0 m (this product requires a trunk cable)**
 - 3 leads snap (BIOECG.RA.S3),
 - 3 leads clamp (BIOECG.RA.C3),
 - 5 leads snap (BIOECG.RA.S5),
 - 5 leads clamp (BIOECG.RA.C5),
 - **Compatible with:**
 - M2000, M3000, M9100,
- Cable extension available on demand,
- 12 months warranty



Moriya	Options	Code
3 lead snap	3 lead snap	BIOECG.BC.S3F
3 lead clamp	3 lead clamp	BIOECG.BC.C3F
5 lead snap	5 lead snap	BIOECG.BC.S5F
5 lead clamp	5 lead clamp	BIOECG.BC.C5F

Moriya	Options	Code
3 lead snap	3 lead snap	BIOECG.BC.S3D
3 lead clamp	3 lead clamp	BIOECG.BC.C3D
5 lead snap	5 lead snap	BIOECG.BC.S5D
5 lead clamp	5 lead clamp	BIOECG.BC.C5D

Moriya	Options	Code
3 lead snap	3 lead snap	BIOECG.BC.3T
5 lead snap	5 lead snap	BIOECG.BC.5T

Moriya	Options	Code
3 lead snap	3 lead snap	BIOECG.RA.S3
3 lead clamp	3 lead clamp	BIOECG.RA.C3
5 lead snap	5 lead snap	BIOECG.RA.S5
5 lead clamp	5 lead clamp	BIOECG.RA.C5