



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Nihon Kohden ECG Cable (New)

- New,
- Made in Uruguay,
- Code: IEC
- Options:
 - **Cable Length: 3,5 m**
 - 3 leads snap (BIOECG.NK.S3F),
 - 3 leads clamp (BIOECG.NK.C3F),
 - 5 leads snap (BIOECG.NK.S5F),
 - 5 leads clamp (BIOECG.NK.C5F),
 - **Compatible with:**
 - Life Scope 6, Life Scope 8, Life Scope 9, Life Scope 11, 6102, 6205, 6120, OEC 6102A (S/N < 20156), OEC-6105 (S/N < 20226), OEC-6120 A/J/K,
 - **Cable Length: 4,0 m**
 - 3 leads snap (BIOECG.NK.S3D),
 - 3 leads clamp (BIOECG.NK.C3D),
 - 5 leads snap (BIOECG.NK.S5D),
 - 5 leads clamp (BIOECG.NK.C5D),
 - **Compatible with:**
 - Life Scope 6, Life Scope 8, Life Scope 9, Life Scope 11, 6102, 6205, 6120, OEC 6102A (S/N < 20156), OEC-6105 (S/N < 20226), OEC-6120 A/J/K,
 - **Cable Length: 3,0 m (this product requires a lead set)**
 - 3 leads (BIOECG.NK.3T),
 - 5 leads (BIOECG.NK.5T),
 - **Compatible with:**
 - Life Scope 6, Life Scope 8, Life Scope 9, Life Scope 11, 6102, 6205, 6120, OEC 6102A (S/N < 20156),OEC-6105 (S/N < 20226), OEC-6120 A/J/K,
 - **Cable Length: 1,0 m (this product requires a trunk cable)**
 - 3 leads snap (BIOECG.RA.S3),
 - 3 leads clamp (BIOECG.RA.C3),
 - 5 leads snap (BIOECG.RA.S5),
 - 5 leads clamp (BIOECG.RA.C5),
 - **Compatible with:**
 - Life Scope 6, Life Scope 8, Life Scope 9, Life Scope 11, 6102, 6205, 6120, OEC 6102A (S/N < 20156), OEC-6105 (S/N < 20226), OEC-6120 A/J/K,
 - **Cable Length: 3,5 m**
 - 3 leads snap (BIOECG.NKV.S3F),
 - 3 leads clamp (BIOECG.NKV.C3F),
 - 5 leads snap (BIOECG.NKV.S5F),
 - 5 leads clamp (BIOECG.NKV.C5F),
 - **Compatible with:**
 - Merlin 78352C, Merlin 78354C, Merlin 78834C, M1001A, M1002A, Viridia, Omnicare CMA, ML65A, 66A, Heartstream and Codemaster XL, Tec-7621K, Tec-7631K, Tec-7721C, Tec-7731C, Tec-7700, Tec-5521, Tec-5531,
 - **Cable Length: 4,0 m**
 - 3 leads snap (BIOECG.NKV.S3D),
 - 3 leads clamp (BIOECG.NKV.C3D),
 - 5 leads snap (BIOECG.NKV.S5D),
 - 5 leads clamp (BIOECG.NKV.C5D),
 - **Compatible with:**
 - Merlin 78352C, Merlin 78354C, Merlin 78834C, M1001A, M1002A, Viridia, Omnicare CMA, ML65A, 66A, Heartstream and Codemaster XL, Tec-7621K, Tec-7631K, Tec-7721C, Tec-7731C, Tec-7700, Tec-5521, Tec-5531,
 - **Cable Length: 3,0 m (this product requires a lead set)**
 - 3 leads (BIOECG.NKV.3T),
 - 5 leads (BIOECG.NKV.5T),
 - **Compatible with:**
 - Merlin 78352C, Merlin 78354C, Merlin 78834C, M1001A, M1002A, Viridia, Omnicare CMA, ML65A, 66A, Heartstream and Codemaster XL, Tec-7621K, Tec-7631K, Tec-7721C, Tec-7731C, Tec-7700, Tec-5521, Tec-5531,
 - **Cable Length: 1,0 m (this product requires a trunk cable)**
 - 3 leads snap (BIOECG.RA.S3),
 - 3 leads clamp (BIOECG.RA.C3),
 - 5 leads snap (BIOECG.RA.S5),
 - 5 leads clamp (BIOECG.RA.C5),
 - **Compatible with:**
 - Merlin 78352C, Merlin 78354C, Merlin 78834C, M1001A, M1002A, Viridia, Omnicare CMA, ML65A, 66A, Heartstream and Codemaster XL, Tec-7621K, Tec-7631K, Tec-7721C, Tec-7731C, Tec-7700, Tec-5521, Tec-5531,
 - **Cable Length: 3,5 m**
 - 3 leads snap (BIOECG.NKT.S3F),
 - 3 leads clamp (BIOECG.NKT.C3F),
 - 5 leads snap (BIOECG.NKT.S5F),
 - 5 leads clamp (BIOECG.NKT.C5F),
 - **Compatible with:**
 - BSM-2300 Life Scope I, BSM-2301A, BSM-2301K, BSM-2351A, BSM-4100, BSM-4102, BSM-4102A, Life Scope P, BSM-5100 Life Scope A, Life Scope TR, BSM-6000 series, BSM-6301K, BSM-6501K, BSM-6701K, Life Scope J, BSM-9101, PVM-2701, BSM-9510 Life Scope M, BSM-9500 Life Scope S.,
 - **Cable Length: 4,0 m**
 - 3 leads snap (BIOECG.NKT.S3D),
 - 3 leads clamp (BIOECG.NKT.C3D),
 - 5 leads snap (BIOECG.NKT.S5D),
 - 5 leads clamp (BIOECG.NKT.C5D),
 - **Compatible with:**
 - BSM-2300 Life Scope I, BSM-2301A, BSM-2301K, BSM-2351A, BSM-4100, BSM-4102, BSM-4102A, Life Scope P, BSM-5100 Life Scope A, Life Scope TR, BSM-6000 series, BSM-6301K, BSM-6501K, BSM-6701K, Life Scope J, BSM-9101, PVM-2701, BSM-9510 Life Scope M, BSM-9500 Life Scope S,
 - **Cable Length: 3,0 m (this product requires a lead set)**
 - 3 leads (BIOECG.NKT.3T),
 - 5 leads (BIOECG.NKT.5T),
 - **Compatible with:**
 - BSM-2300 Life Scope I, BSM-2301A, BSM-2301K, BSM-2351A, BSM-4100, BSM-4102, BSM-4102A, Life Scope P, BSM-5100 Life Scope A, Life Scope TR, BSM-6000 series, BSM-6301K, BSM-6501K, BSM-6701K, Life Scope J, BSM-9101, PVM-2701, BSM-9510 Life Scope M, BSM-9500 Life Scope S,
 - **Cable Length: 1,0 m (this product requires a trunk cable)**
 - 3 leads snap (BIOECG.RA.S3),
 - 3 leads clamp (BIOECG.RA.C3),
 - 5 leads snap (BIOECG.RA.S5),
 - 5 leads clamp (BIOECG.RA.C5),
 - **Compatible with:**
 - BSM-2300 Life Scope I, BSM-2301A, BSM-2301K, BSM-2351A, BSM-4100, BSM-4102, BSM-4102A, Life Scope P, BSM-5100 Life Scope A, Life Scope TR, BSM-6000 series, BSM-6301K, BSM-6501K, BSM-6701K, Life Scope J, BSM-9101, PVM-2701, BSM-9510 Life Scope M, BSM-9500 Life Scope S,
 - **Cable Length: 3,0 m**
 - 10 lead banana 4 mm (BIOECG.NK.P10F),
 - **Compatible with:**
 - Cardiofax 651, 6581, 8110, 9010, 9110, 9620, Cardiofax Q, 1350, ECG-1250A, ECG-9620, ECG-9010, ECG-9020, ECG-9022, ECG-9110, ECG-9130, ECG-9132, ECG-9110K, ECG-9130K, Cardioline (new line), Del Mar Reynolds Voyager 12; Newtech 1206; Carewell 1103, Carewell AR 1200ADV,
 - **Cable Length: 3,0 m**
 - 10 lead banana 4 mm (BIOECG.NK2.P10F),
 - **Compatible with:**
 - Cardiofax/Q-Cardiofax, 6151, 6353, 6851, 8110, 8270, 8330, 8330A, 8350, 8370, 8420, 8423, 8450, 8453, 8820, 8823, 8830, 9020, 9022, 9130, 9132, Escape 12, Edan Smart ECG, Danatech, Dr. Lee 120B
 - **Cable extension available on demand,**
 - 12 months warranty



Włose Kubiki

Opis	Waga	Wartość
Włose Kubiki	1000	1000
Włose Kubiki	1000	1000
Włose Kubiki	1000	1000
Włose Kubiki	1000	1000

Włose Kubiki

Opis	Waga	Wartość
Włose Kubiki	1000	1000
Włose Kubiki	1000	1000
Włose Kubiki	1000	1000
Włose Kubiki	1000	1000

Włose Kubiki

Opis	Waga	Wartość
Włose Kubiki	1000	1000
Włose Kubiki	1000	1000
Włose Kubiki	1000	1000
Włose Kubiki	1000	1000

Włose Kubiki

Opis	Waga	Wartość
Włose Kubiki	1000	1000
Włose Kubiki	1000	1000
Włose Kubiki	1000	1000
Włose Kubiki	1000	1000

Włose Kubiki

Opis	Waga	Wartość
Włose Kubiki	1000	1000
Włose Kubiki	1000	1000
Włose Kubiki	1000	1000
Włose Kubiki	1000	1000

Włose Kubiki

Opis	Waga	Wartość
Włose Kubiki	1000	1000
Włose Kubiki	1000	1000
Włose Kubiki	1000	1000
Włose Kubiki	1000	1000

Włose Kubiki

Opis	Waga	Wartość
Włose Kubiki	1000	1000
Włose Kubiki	1000	1000
Włose Kubiki	1000	1000
Włose Kubiki	1000	1000

Włose Kubiki

Opis	Waga	Wartość
Włose Kubiki	1000	1000
Włose Kubiki	1000	1000
Włose Kubiki	1000	1000
Włose Kubiki	1000	1000

Włose Kubiki

Opis	Waga	Wartość
Włose Kubiki	1000	1000
Włose Kubiki	1000	1000
Włose Kubiki	1000	1000
Włose Kubiki	1000	1000