



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Progetti ECG Cable (New)

- New,
- Made in Uruguay,
- 10 lead cable,
- Code: IEC
- Plug: banana 4 mm,
- **Cable Length 3,0 m** (cable extension available on demand),
- Compatible with:
 - EPG Smart, EPG Pocket, EPG Junior, EPG View, EPG 6 View, EPG Plus, EPG 612, EPG 13,
- 12 months warranty

