



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Schiller ECG Cable (New)

- New,
- Made in Uruguay,
- Code: IEC
- Options:
 - **Cable Length: 3,0 m**
 - 10 lead banana (BIOECG.SC.P10F),
 - **Compatible with:**
 - AT3, AT6, CS6, AT5, AT10, AT60, CS-200; Bionet Cardiotouch,
 - **Cable Length: 3,0 m**
 - 10-odpr. banan (BIOECG.SC2.P10F),
 - **Compatible with:**
 - AT-1, AT-2, AT-2 Light, AT-2 Plus, AT-4, AT-10, AT-10 Plus, AT-60, AT-101, AT-102, AT-102-Plus, AT-104, AT-110, CS-100, CS-200, MS-2010, MS-2015; Bionet: CardioCare 2000; Esaote: P80 Series 40893, Biosound MyLab 30 CD, 30 CV; Cardiette: AR1200 adv; Welch Allyn CP10, CP20,
- **Cable extension available on demand,**
- 12 months warranty

