



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Spacelabs ECG Cable (New)

- New,
- Made in Uruguay,
- Code: IEC
- Options:
 - **Cable Length: 3,5 m**
 - 3 leads snap (BIOECG.BC.S3F),
 - 3 leads clamp (BIOECG.BC.C3F),
 - 5 leads snap (BIOECG.BC.S5F),
 - 5 leads clamp (BIOECG.BC.C5F),
 - **Compatible with:**
 - 500 / 600 / 900 Series (w/o Respiration), 90300 / 90400 / 90600 / 90700 Series, Tek-208 / 408 / 511 / 512 514 / 521, PC Bedsi, PCMS,
 - **Cable Length: 4,0 m**
 - 3 leads snap (BIOECG.BC.S3D),
 - 3 leads clamp (BIOECG.BC.C3D),
 - 5 leads snap (BIOECG.BC.S5D),
 - 5 leads clamp (BIOECG.BC.C5D),
 - **Compatible with:**
 - 500 / 600 / 900 Series (w/o Respiration), 90300 / 90400 / 90600 / 90700 Series, Tek-208 / 408 / 511 / 512 514 / 521, PC Bedsi, PCMS,
 - **Cable Length: 3,0 m (this product requires a lead set)**
 - 3 leads (BIOECG.BC.3T),
 - 5 leads (BIOECG.BC.5T),
 - **Compatible with:**
 - 500 / 600 / 900 Series (w/o Respiration), 90300 / 90400 / 90600 / 90700 Series, Tek-208 / 408 / 511 / 512 514 / 521, PC Bedsi, PCMS,
 - **Cable Length: 1,0 m (this product requires a trunk cable)**
 - 3 leads snap (BIOECG.RA.S3),
 - 3 leads clamp (BIOECG.RA.C3),
 - 5 leads snap (BIOECG.RA.S5),
 - 5 leads clamp (BIOECG.RA.C5),
 - **Compatible with:**
 - 500 / 600 / 900 Series (w/o Respiration), 90300 / 90400 / 90600 / 90700 Series, Tek-208 / 408 / 511 / 512 514 / 521, PC Bedsi, PCMS,
 - **Cable Length: 3,5 m**
 - 3 leads snap (BIOECG.SP.S3F),
 - 3 leads clamp (BIOECG.SP.C3F),
 - 5 leads snap (BIOECG.SP.S3F),
 - 5 leads clamp (BIOECG.SP.C3F),
 - **Compatible with:**
 - Ultraview 90309, 90367, 90492, 90496, 91220, 91496,
 - **Cable Length: 4,0 m**
 - 3 leads snap (BIOECG.SP.S3D),
 - 3 leads clamp (BIOECG.SP.C3D),
 - 5 leads snap (BIOECG.SP.S5D),
 - 5 leads clamp (BIOECG.SP.C5D),
 - **Compatible with:**
 - Ultraview 90309, 90367, 90492, 90496, 91220, 91496,
 - **Cable Length: 3,0 m (this product requires a lead set)**
 - Kabel główny 3- odpr. (BIOECG.SP.3T),
 - Kabel główny 5-odpr. (BIOECG.SP.5T),
 - **Compatible with:**
 - Ultraview 90309, 90367, 90492, 90496, 91220, 91496,
 - **Cable Length: 1,0 m (this product requires a trunk cable)**
 - 3 leads snap (BIOECG.RA.S3),
 - 3 leads clamp (BIOECG.RA.C3),
 - 5 leads snap (BIOECG.RA.S5),
 - 5 leads clamp (BIOECG.RA.C5),
 - **Compatible with:**
 - Ultraview 90309, 90367, 90492, 90496, 91220, 91496,
- **Cable extension available on demand,**
- 12 months warranty

[illegible][illegible]