



**PRAISTON Sp. z o.o.**  
ul. Górska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## TEB ECG Cable (New)

- New,
- Made in Uruguay,
- Code: IEC
- Options:
  - **Cable Length: 3,5 m**
    - 3 leads snap (BIOECG.TE.S3F),
    - 3 leads clamp (BIOECG.TE.C3F),
    - 5 leads snap (BIOECG.TE.S5F),
    - 5 leads clamp (BIOECG.TE.C5F),
  - **Compatible with:**
    - M10, M10B, Med1, M10A, M10C,
  - **Cable Length: 4,0 m**
    - 3 leads snap (BIOECG.TE.S3D),
    - 3 leads clamp (BIOECG.TE.C3D),
    - 5 leads snap (BIOECG.TE.S5D),
    - 5 leads clamp (BIOECG.TE.C5D),
  - **Compatible with:**
    - M10, M10B, Med1, M10A, M10C,
  - **Cable Length: 3,0 m (this product requires a lead set)**
    - 3 leads (BIOECG.TE.3T),
    - 5 leads (BIOECG.TE.5T)
  - **Compatible with :**
    - M10, M10B, Med1, M10A, M10C,
  - **Cable Length: 1,0 m (this product requires a trunk cable)**
    - 3 leads snap (BIOECG.RA.S3),
    - 3 leads clamp (BIOECG.RA.C3),
    - 5 leads snap (BIOECG.RA.S5),
    - 5 leads clamp (BIOECG.RA.C5),
  - **Compatible with:**
    - M10, M10B, Med1, M10A, M10C
- **Cable extension available on demand,**
- 12 months warranty



| TEB                           | Options       | Code   |
|-------------------------------|---------------|--------|
| 3,5 m cable with 3 or 5 leads | 3 leads snap  | BIOS3F |
|                               | 3 leads clamp | BIOC3F |
|                               | 5 leads snap  | BIOS5F |
|                               | 5 leads clamp | BIOC5F |

| TEB                           | Options       | Code   |
|-------------------------------|---------------|--------|
| 4,0 m cable with 3 or 5 leads | 3 leads snap  | BIOS3D |
|                               | 3 leads clamp | BIOC3D |
|                               | 5 leads snap  | BIOS5D |
|                               | 5 leads clamp | BIOC5D |

| TEB                           | Options | Code   |
|-------------------------------|---------|--------|
| 3,0 m cable with 3 or 5 leads | 3 leads | BIOS3T |
|                               | 5 leads | BIOS5T |

| TEB                           | Options       | Code   |
|-------------------------------|---------------|--------|
| 1,0 m cable with 3 or 5 leads | 3 leads snap  | BIOS3S |
|                               | 3 leads clamp | BIOC3S |
|                               | 5 leads snap  | BIOS5S |
|                               | 5 leads clamp | BIOC5S |