



**PRAISTON Sp. z o.o.**  
 ul. Górowska 32  
 64-100 Leszno  
 Polska

Biuro: +48 65 527 01 67  
 Serwis: +48 782 844 000  
 Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
 Krzysztof Wybieralski  
 +48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## Welch Allyn ECG Cable (New)

- New,
- Made in Uruguay,
- Code: IEC
- Options:
  - **Cable Length: 3,5 m**
    - 3 leads snap (BIOECG.BC.S3F),
    - 3 leads clamp (BIOECG.BC.C3F),
    - 5 leads snap (BIOECG.BC.S5F),
    - 5 leads clamp (BIOECG.BC.C5F),
  - **Compatible with:**
    - Propaq 100/200 Series, Atlas Propaq Encore, Propaq CS, Ultra Smart Protocol,
  - **Cable length: 4,0 m**
    - 3 leads snap (BIOECG.BC.S3D),
    - 3 leads clamp (BIOECG.BC.C3D),
    - 5 leads snap (BIOECG.BC.S5D),
    - 5 leads clamp (BIOECG.BC.C5D),
  - **Compatible with:**
    - Propaq 100/200 Series, Atlas Propaq Encore, Propaq CS, Ultra Smart Protocol,
  - **Cable Length: 3,0 m (this product requires a lead set)**
    - 3 leads (BIOECG.BC.3T),
    - 5 leads (BIOECG.BC.5T),
  - **Compatible with:**
    - Propaq 100/200 Series, Atlas Propaq Encore, Propaq CS, Ultra Smart Protocol,
  - **Cable Length: 1,0 m (this product requires a trunk cable)**
    - 3 leads snap (BIOECG.RA.S3),
    - 3 leads clamp (BIOECG.RA.C3),
    - 5 leads snap (BIOECG.RA.S5),
    - 5 leads clamp (BIOECG.RA.C5),
  - **Compatible with:**
    - Propaq 100/200 Series, Atlas Propaq Encore, Propaq CS, Ultra Smart Protocol,
- **Cable extension available on demand,**
- 12 months warranty

