



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Zoll ECG cable (New)

- New,
- Made in Uruguay,
- Code: IEC
- Options:
 - **Cable Length: 3,5 m**
 - 3 leads snap (BIOECG.BC.S3F),
 - 3 leads clamp (BIOECG.BC.C3F),
 - 5 leads snap (BIOECG.BC.S5F),
 - 5 leads clamp (BIOECG.BC.C5F),
 - **Compatible with:**
 - D900, PD1200, PD1400, PD1600, PD1700, PD2000,
 - **Cable Length: 4,0 m**
 - 3 leads snap (BIOECG.BC.S3D),
 - 3 leads clamp (BIOECG.BC.C3D),
 - 5 leads snap (BIOECG.BC.S5D),
 - 5 leads clamp (BIOECG.BC.C5D),
 - **Compatible with:**
 - D900, PD1200, PD1400, PD1600, PD1700, PD2000,
 - **Cable Length: 3,0 m (this product requires a lead set)**
 - 3 leads (BIOECG.BC.3T),
 - 5 leads (BIOECG.BC.5T),
 - **Compatible with:**
 - D900, PD1200, PD1400, PD1600, PD1700, PD2000,
 - **Cable Length: 1,0 m (this product requires a trunk cable)**
 - 3 leads snap (BIOECG.RA.S3),
 - 3 leads clamp (BIOECG.RA.C3),
 - 5 leads snap (BIOECG.RA.S5),
 - 5 leads clamp (BIOECG.RA.C5),
 - **Compatible with:**
 - D900, PD1200, PD1400, PD1600, PD1700, PD2000
- **Cable extension available on demand,**
- 12 months warranty



Zoll	System	Code
ECG lead cable 3,5 m 3 leads	3 leads snap	BIOECG.BC.S3F
ECG lead cable 3,5 m 3 leads	3 leads clamp	BIOECG.BC.C3F
ECG lead cable 3,5 m 5 leads	5 leads snap	BIOECG.BC.S5F
ECG lead cable 3,5 m 5 leads	5 leads clamp	BIOECG.BC.C5F

Zoll	System	Code
ECG lead cable 4,0 m 3 leads	3 leads snap	BIOECG.BC.S3D
ECG lead cable 4,0 m 3 leads	3 leads clamp	BIOECG.BC.C3D
ECG lead cable 4,0 m 5 leads	5 leads snap	BIOECG.BC.S5D
ECG lead cable 4,0 m 5 leads	5 leads clamp	BIOECG.BC.C5D

Zoll	System	Code
ECG lead cable 3,0 m 3 leads	3 leads	BIOECG.BC.3T
ECG lead cable 3,0 m 5 leads	5 leads	BIOECG.BC.5T

Zoll	System	Code
ECG lead cable 1,0 m 3 leads	3 leads snap	BIOECG.RA.S3
ECG lead cable 1,0 m 3 leads	3 leads clamp	BIOECG.RA.C3
ECG lead cable 1,0 m 5 leads	5 leads snap	BIOECG.RA.S5
ECG lead cable 1,0 m 5 leads	5 leads clamp	BIOECG.RA.C5