



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Gynaecological Patient Bed TEYCO MED FRANCY 4 (New)

- New,
- Made in Italy,
- Equipped with 4 powerful actuators,
- Made entirely of steel painted with epoxy,
- Electrical Movements:
 - Height adjustment,
 - Leg Plates,
 - Back Rest,
 - Trendelenburg,
- On casters (with break),
- Rotation of the seat to the left and to the right at 45° degrees with manual locking,
- Technical Specifications:
 - Power input: 230 V – 50 Hz,
 - Power: 24 V,
 - Height: 10000 N,
 - Back Rest/Foot Rest: 6000 N
 - Trendelenburg: 2500 N,
 - Maximum load: 220 - 230 kg,
 - Maximum Power: 200 VA,
 - IP66 protection standard,
 - Weight: 150 kg,
 - Sitting height max: 910 mm,
 - Sitting height min: 518 mm,
- **Additional Equipment (extra paid):**
 - Headrest,
 - Paper dispenser,
 - IV Pole,
 - Side rail for additional equipment,
 - Accessory tray,
 - Upholstery (blue color),
 - Two-colour upholstery (Back/Sitting + Foot Rest),
 - Tri-colour upholstery (Back + Sitting + Foot Rest),
 - Foot switch (Up /Down),
 - 4 gummy wheels with metal frame + breaks,
- Technical Passport,
- 24-month warranty

