



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Gynaecological Patient Bed TEYCO MED ELLY03 (New)

- New,
- Made in Italy,
- Equipped with 3 powerful actuators,
- Made entirely of steel painted with epoxy,
- Electrical Movements:
 - Height adjustment,
 - Trendelenburg,
 - Back Rest,
- Leg Plates regulated by gas spring,
- Stainless Steel Irrigation bowl,
- Technical Specifications:
 - Power input: 220/230 V – 50 Hz,
 - Power: 24 V,
 - Height: 6000 N,
 - Back Rest: 3000 N,
 - Maximum Load: 180 kg,
 - Maximum Power: 120 VA,
 - Fuse: 1,6 + 6 A,
 - IP44 protection standard,
 - Weight: 90 kg,
 - Base Length: 800 mm,
 - Sitting height max: 840 mm,
 - Sitting height min: 480 mm,
- **Additional Equipment (extra paid):**
 - Telescopic pillar (additional or exchanged for leg plates),
 - Colposcopy support,
 - Headrest,
 - Paper dispenser,
 - Electrical Back Rest,
 - Electrical Leg Plates,
 - IV Pole,
 - Knee Clutch,
 - Casters with break system,
 - Knee Clutch cover,
 - Metal telescopic pillar model RING,
 - Chrome Base,
 - Transparent Glass Base: 19 mm,
- Technical Passport,
- 24-month warranty

