



Examination Chair TEYCO MED EXPEDIA EB01 (Expedia basic) (New)

- New,
- Made in Italy,
- Equipped with 1 powerful actuators,
- Folded footstool,
- Sitting: adjustable angle 45 ° left/right,
- Electrical Movements: Height,
- Back rest and footstool regulated by gas spring,
- Technical Specifications:
 - Power Input: 230 V – 50 Hz,
 - Power: 24 V,
 - Height: 6000 N,
 - Gas spring: 800 N,
 - Max Load: 180 kg,
 - Max power: 120 VA,
 - Fuse: 1,6 + 6 A,
 - IP44 protection standard,
 - Weight: 75 kg,
 - Base length: 580 x 740 mm,
 - Sitting height max: 980 mm,
 - Sitting height min: 650 mm,
- **Additional Equipment (extra paid):**
 - Headrest,
 - IV Pole,
 - Casters with breaks, metal frame,
- Technical Passport,
- 24-month warranty

