



## Phlebotomy Chair TEYCO MED LORY Y01 (New)

- New
- Made in Italy
- Equipped with 2 powerful actuators,
- Electrical Movements:
  - Height adjustment,
  - Back Rest,
- Technical Data:
  - Power supply: 230V – 50 Hz,
  - Power: 24 V,
  - Back Rest/Foot Rest: 6000 N,
  - Gas Spring: 800 N,
  - Max Load: 180 kg,
  - Max Power: 120 VA,
  - Fuse: 1,6 + 6 A,
  - IP44 protection standard,
  - Weight: 75 kg,
  - Base Length: 580 x 740 mm,
  - Sitting height max: 980 mm,
  - Sitting height min: 650 mm,
- **Additional Equipment (paid):**
  - Headrest,
  - Paper dispenser,
  - IV Pole,
  - 4 gummy wheels with metal frame + breaks,
- Technical Passport,
- 24 months warranty

