



Examination Chair TEYCO MED FRANCY B3 (New)

- New,
- Made in Italy,
- Equipped with 3 powerful actuators,
- Made entirely of steel painted with epoxy,
- Electrical Movements:
 - Height adjustment,
 - Leg Plates
 - Back Rest,
 - Trendelenburg,
- On casters (with break),
- Technical Specifications:
 - Power input: 230 V – 50 Hz,
 - Power: 24 V,
 - Height: 10000 N,
 - Back Rest/Foot Rest: 6000 N
 - Trendelenburg: 2500 N,
 - Maximum load: 220-230 kg,
 - Maximum Power: 200 VA,
 - IP66 protection standard,
 - Weight: 150 kg,
 - Sitting height max: 910 mm,
 - Sitting height min: 518 mm,
- **Additional Equipment (extra paid):**
 - Headrest,
 - Paper dispenser,
 - IV Pole,
 - Accessory tray,
 - Actuator,
 - Side rail for additional equipment,
 - Upholstery (blue color),
 - Two-colour upholstery (Back/Sitting + Foot Rest),
 - Tri-colour upholstery (Back + Sitting + Foot Rest),
 - Foot switch (Up /Down),
- Technical Passport,
- 24-month warranty

