



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **Examination Table TEYCO MED FRANCY O4 (New)**

- New,
- Made in Italy,
- Equipped with 4 powerful actuators,
- Made entirely of steel painted with epoxy,
- Electrical Movements:
  - Height adjustment,
  - Leg Plates,
  - Back Rest,
  - Trendelenburg,
- On casters (with break),
- Rotation of the seat to the left and to the right at 45° degrees with manual locking,
- Technical Specifications:
  - Power input: 230 V – 50 Hz,
  - Power: 24 V,
  - Height: 10000 N,
  - Back Rest/Foot Rest: 6000 N
  - Trendelenburg: 2500 N,
  - Maximum load: 220 - 230 kg,
  - Maximum Power: 200 VA,
  - IP66 protection standard,
  - Weight: 150 kg,
  - Sitting height max: 910 mm,
  - Sitting height min: 518 mm,
- **Additional Equipment (extra paid):**
  - Headrest,
  - Paper dispenser,
  - IV Pole,
  - Side rail for additional equipment,
  - Accessory tray,
  - Upholstery (blue color),
  - Two-colour upholstery (Back/Sitting + Foot Rest),
  - Tri-colour upholstery (Back + Sitting + Foot Rest),
  - Foot switch (Up /Down),
  - 4 gummy wheels with metal frame + breaks,
- Technical Passport,
- 24-month warranty,
- Financing:
  - Installments,
  - Leasing,
  - Loan.

