



Examination Chair TEYCO MED IVY1 (New)

- New,
- Made in Italy,
- Equipped with 1 powerful actuators,
- Made entirely of steel painted with epoxy,
- Electrical Movements:
 - Height adjustment,
- Back Rest and Leg Plates regulated by gas spring,
- Stainless Steel Irrigation bowl,
- Technical Specifications:
 - Power input: 230 V – 50 Hz,
 - Power: 24 V,
 - Height: 6000 N,
 - Gas Spring: 800 N,
 - Maximum Load: 180 kg,
 - Maximum Power: 120 VA,
 - Fuse: 1,6 + 6 A,
 - IP44 protection standard,
 - Weight: 75 kg,
 - Base Length: 580 x 740 mm,
 - Sitting height max: 950 mm,
 - Sitting height min: 500 mm,
- **Additional Equipment (extra paid):**
 - Telescopic pillar (additional or exchanged for leg plates),
 - Colposcopy support,
 - Electrical Back Rest,
 - Electrical Leg Plates,
 - Headrest,
 - Paper dispenser,
 - IV Pole,
 - Knee Clutch,
 - Casters with break system,
 - Knee Clutch cover,
 - Metal telescopic pillar model RING,
- Technical Passport,
- 24-month warranty

