



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Internal Electrodes for MEDTRONIC (PHYSIO-CONTROL) LIFEPAK 9 (Reconditioned)

- Condition: Reconditioned - Very good,
- Made in USA,
- Open heart electrodes,
- Cable length: 300 cm,
- Electrode diameter: 50 mm,
- Technical Passport,
- Warranty: 30-days

