



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



BLOODY 03 blood collection chair (New)

New,
Made in Poland,
The design makes it easy to carry out the injection procedure on the left or right forearm of the person subject to the procedure,
There is a free option to order upholstery in a colour of your choice, except for "PRESTIGE" colours,
The chair has a stable steel construction with powder coating, adjustable backrest and two armrests adjustable in all planes. Backrest and footrest are adjusted by means of an electric actuator.
It is used in doctor's surgeries, health clinics, hospitals and laboratories,
The armchair can be fitted for an additional charge:

Ordinary headrest,
Laryngological headrest,
Profiled upholstery,
Technical parameters:
Height 101 cm,
Width 50 cm, total width 71 cm,
Depth 81 cm,
Overall length when unfolded with footrest 170 cm,
Seat height 50 cm,
own weight 35 kg,
Maximum load 150 kg,
Headrest height adjustment 15 cm
Height adjustment of armrests 25 cm,
Backrest angle adjustment 0°- 85°,
Headrest angle adjustment 0°- 85°,
It has a current Technical Passport,

Warranty: 24 months,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN NEDERLANDS

МЫ ГОВОРИМ ПО РУССКИ

МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ

WE SHIP WORLDWIDE

