



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Waiting Room Bench (L2TPS1UT) 2 seater (New)

- New,
- Made in Poland,
- Made from high quality plywood,
- Iron frames painted with powder coating,
- With automatic hinges
- Floor-mounted
- Seat with upholstery
- Dimensions:
 - Width: 117 cm,
 - Depth: 50,5 cm,
 - Height: 80 cm ,
- Available color configurations:
 - Plywood:
 - grey + natural fringes
 - graphite + natural fringes
 - white + natural fringes
 - natural oak
 - Frame:
 - grey,
 - graphite,
 - white
- 24 months warranty

