



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Rehabilitation Bed DREAM 325 (Linak) (New)

New

Made in Poland

The bed has a four-section sleep surface

Backrest and thigh sections auto-contouring function

Sleep surface height adjustment from 40cm to 80cm

All bed functions are controlled electrically with a remote control

Option to block particular positions of the bed (Linak)

Powder-varnished metal elements

Wooden elements in beech colour

Four antistatic casters with blocking option

Extension arm and side rails included as standard

Bed fitted with LINAK® actuators

Technical parameters

LINAK drive and control

Safety working load: 200 kg

Maximum user weight: 165 kg

Total weight: 87.0 kg

External dimensions: length 218 cm, width 102 cm

Backrest section adjustment: 0 ÷ 70°

Adjustment of lower backrest section: 0 ÷ 20

Sleep surface height adjustment with electric motors: 40 - 80 cm,

Clearance under the bed: approx. 24 cm,

Castors: 100mm with a brake,

The bed has a valid technical passport

Warranty: 24 months,

Financing options (Poland only): instalments, leasing, loan.

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN NEDERLANDS

МЫ ГОВОРИМ ПО РУССКИ

МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ

WE SHIP WORLDWIDE

