



Medicare Standard 2 Section (SM2570) (New)

- New,
- Made in UK,
- Electric Elevation, electric backrest,
- Retract Wheel Base,
- Remote Control,
- Vertical lifting design with stability throughout,
- Anti-bacterial powder coated framework,
- TiMotion electric linear actuators,
- Technical Specifications:
 - Length: 194 cm,
 - Width: 63 cm,
 - Safe working load and lifting capacity of 240 kg (528 lbs),
 - Adjustable angle electric backrest from -25° to 75°,
 - Excellent height range from 45 cm to 98 cm,
- **Additional Equipment (extra paid):**
 - Perimeter foot switch operation for electric models (6001),
 - Side support rails (Pair) (6002),
 - Side support rails cushions (Pair) (6002P),
 - Breathing hole and plug for prone treatments (6004),
 - Paper Roll Handle (6005),
 - Removable head/knee support (6006),
 - Removable headrest - pillow (6007),
 - Protective Base Shrouding (6053),
 - IV Pole (6054),
 - Battery Back-up (6057),
 - Foot switch (6200),
- Technical Passport,
- 24-month warranty for electric

