



## Medicare Mobile Treatment 2 (SM2660) (New)

- New,
- Made in UK,
- Fully electric height adjustment,
- Adjustable angle of the backrest controlled manually by a gas spring,
- Large Wheel Mobile Base Design for easy movement,
- Compact frame design,
- Basic equipment:
  - Paper roll handle (6005),
  - Fitted fold up/down Side support rails (6002),
- Anti-bacterial powder coated framework,
- TiMotion electric linear actuators,
- Technical Specifications:

- Length: 194 cm,
- Width: 63 cm made from foam,
- Safe working load and lifting capacity of 240 kg (528 lbs),
- Adjustable angle electric backrest from -25° to 75°,
- Excellent height range from 46 cm to 99 cm,
- Wheels Tente® Ø125 mm,

### • Additional Equipment (extra paid):

- Cushions for Side Rails – a pair (6002P),
- Breathing hole and plug (6004),
- Removable head/knee support (6006),
- Removable headrest-pillow (6007),
- Base cover (6053),
- IV Pole (6054),
- Foot switch (6200),
- Technical Passport,
- 24-month warranty for electric

