



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Medicare Mobile Treatment 2 (SM2670) (New)

- New,
- Made in UK,
- Fully electric height adjustment as well as backrest,
- Large Wheel Mobile Base Design for easy movement,
- Compact frame design,
- Basic equipment:
 - Paper roll handle (6005),
 - Fitted fold up/down Side support rails (6002),
- Anti-bacterial powder coated framework,
- TiMotion electric linear actuators,
- Technical Specifications:
 - Length: 194 cm,
 - Width: 63 cm made from foam,
 - Safe working load and lifting capacity of 240 kg (528 lbs),
 - Adjustable angle electric backrest from -25° to 75°,
 - Excellent height range from 46 cm to 99 cm,
 - Wheels Tente® Ø125 mm,
- **Additional Equipment (extra paid):**
 - Cushions for Side Rails – a pair (6002P),
 - Breathing hole and plug (6004),
 - Removable head/knee support (6006),
 - Removable headrest-pillow (6007),
 - Base cover (6053),
 - IV Pole (6054),
 - Double Foot switch – (6207),
- Technical Passport,
- 24-month warranty for electric

