




Medicare Mobile Treatment 3 (SM3660) (New)

- New,
- Made in UK,
- Fully electric height adjustment,
- Adjustable angle of footrest and backrest controlled by a gas spring,
- Large Wheel Mobile Base Design for easy movement,
- Compact frame design,
- Basic equipment:
 - Paper roll handle (6005),
 - Side support rails - a pair (6003),
- Anti-bacterial powder coated framework,
- TiMotion electric linear actuators,
- Technical Specifications:
 - Length: 191 cm,
 - Width: 63 cm made from foam,
 - Safe working load and lifting capacity of 240 kg (528 lbs),
 - Adjustable angle electric backrest from 0° to 85°,
 - Adjustable angle footrest from -25° to 75°,
 - Excellent height range from 46 cm to 99 cm,
 - Wheels Tente® Ø125 mm,
- **Additional Equipment (extra paid):**
 - Side support cushions - a pair (6003P),
 - Breathing hole and plug (6004),
 - Removable head/knee support (6006),
 - Removable headrest-pillow (6007),
 - Base cover (6053),
 - IV Pole (6054),
 - Foot switch - (6200),
- Technical Passport,
- 24-month warranty for electric



TAPICERKI STANDARD [W CENIE]



TAPICERKI STANDARD +[ZA DOPŁATK.]

