



PRAISTON Sp. z o.o.
 ul. Górowska 32
 64-100 Leszno
 Polska


Biuro: +48 65 527 01 67
 Serwis: +48 782 844 000
 Fax: +48 65 527 01 67
 E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
 Krzysztof Wybieralski
 +48 882 762 006
kw@praiston.pl



Medicare 3 Section (SM3550) (New)

- New,
- Made in UK,
- Hydraulic Elevation,
- Adjustable angle of footrest controlled by a gas spring,
- Individual locking castors or linked central locking wheel system,
- Compact frame design,
- Anti-bacterial powder coated framework,
- TiMotion electric linear actuators,
- Technical Specifications:
 - Length: 191 cm,
 - Width: 63 cm made from foam,
 - Safe working load and lifting capacity of 240 kg (528 lbs),
 - Adjustable angle electric backrest from 0° to 85°,
 - Adjustable angle footrest from -25° to 75°,
 - Excellent height range from 45 cm to 98 cm,
- **Additional Equipment (extra paid):**
 - Side support rails - a pair (6003),
 - Side support cushions - a pair (6003P),
 - Breathing hole and plug (6004),
 - Paper roll handle (6005),
 - Removable head/knee support (6006),
 - Removable headrest-pillow (6007),
 - Base cover (6053),
 - IV Pole (6054),
- Technical Passport,
- 24-month warranty for electric



TAPICERKI STANDARD [W CENIE]



TAPICERKI STANDARD [ZA DOPŁATĄ]

