



**PRAISTON Sp. z o.o.**  
 ul. Górowska 32  
 64-100 Leszno  
 Polska


Biuro: +48 65 527 01 67  
 Serwis: +48 782 844 000  
 Fax: +48 65 527 01 67  
 E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordynator ds. sprzedaży  
 Krzysztof Wybieralski  
 +48 882 762 006  
[kw@praiston.pl](mailto:kw@praiston.pl)



## Medicare 3 Section (SM3550) (New)

- New,
- Made in UK,
- Hydraulic Elevation,
- Adjustable angle of footrest controlled by a gas spring,
- Individual locking castors or linked central locking wheel system,
- Compact frame design,
- Anti-bacterial powder coated framework,
- TiMotion electric linear actuators,
- Technical Specifications:
  - Length: 191 cm,
  - Width: 63 cm made from foam,
  - Safe working load and lifting capacity of 240 kg (528 lbs),
  - Adjustable angle electric backrest from 0° to 85°,
  - Adjustable angle footrest from -25° to 75°,
  - Excellent height range from 45 cm to 98 cm,
- **Additional Equipment (extra paid):**
  - Side support rails - a pair (6003),
  - Side support cushions - a pair (6003P),
  - Breathing hole and plug (6004),
  - Paper roll handle (6005),
  - Removable head/knee support (6006),
  - Removable headrest-pillow (6007),
  - Base cover (6053),
  - IV Pole (6054),
- Technical Passport,
- 24-month warranty for electric



TAPICERKI STANDARD [W CENIE ]



TAPICERKI STANDARD [ZA DOPŁATĄ]

