



Innovation Standard 3 (MG3670) Section Couch (New)

- New,
- Made in UK,
- Fully electric height adjustment as well as an electric profiling foot and backrest,
- Individual locking castors or linked central locking wheel system,
- Suitable for a wide range of general examination and minor treatment procedures,
- Acrylic capped vac formings to facilitate easy cleaning,
- Anti-bacterial powder coated framework,
- TiMotion electric linear actuators,
- Technical Specifications:

- Length: 200 cm,
- Width: 65 cm made from 40 mm foam,
- Safe working load and lifting capacity of 240 kg (528 lbs),
- Adjustable angle electric backrest from 0° to 80°,
- Adjustable angle footrest from -35° to 10°,
- Excellent height range from 45 cm to 95 cm,
- Wheels - Tente® Ø125 mm,

• **Additional Equipment (extra paid):**

- Side support rails - a pair (6003),
- Side support cushions - a pair (6003P),
- Breathing hole and plug (6004),
- Paper roll handle (6005),
- Removable head/knee support (6006),
- Removable headrest-pillow (6007),
- IV Pole (6054),
- Battery back-up (6057),
- Technical Passport,
- 24-month warranty for electric

