



Examination Table Medicare Multi-Couch - Single Footrest (SM9566) (New)

- New,
- Made in UK,
- Electric Elevation,
- Retractable wheel design with adjustable foot for uneven floors,
- Purpose made precision bearings on all pivot points,
- Clearance beneath frame for mobile hoist,
- Hand switch fitted as standard on electric models,
- Foot sections incorporate anti-trap design,
- TiMotion powerful actuators,
- Technical Data:
 - Length: 180 cm,
 - Width: 65 cm made from foam,
 - Safe working load and lifting capacity (from minimum height) of 240 kg (528 lbs = 37 stone),
 - Adjustable back rest with positive angle of 0 to 80 degrees,
 - Adjustable foot section with angle adjustment of 0 to -85 degrees,
 - Excellent height range from 50 - 105 cm,
- **Additional Equipment (extra paid):**
 - Breathing hole and plug for prone treatments (6004),
 - Paper towel holder (6005),
 - Upper limb support (6006),
 - Removable headrest - pillow (6007),
 - Gyneacology stirrups (6044),
 - Gyneacology Leg support (6045),
 - Short base cover (6052),
 - IV Stand (6054),
 - Battery back-up (6057),
 - Single foot switch (6200),
- Technical Passport,
- 24-month warranty for electric, gas spring

