



Medicare Multi-Couch - Single Footrest (SM9576) (New)

- New,
- Made in UK,
- Electric Elevation,
- Retract Wheel Base,
- Vertical lifting design with stability throughout,
- Anti-bacterial powder coated framework,
- TiMotion electric linear actuators,
- Technical Specifications:
 - Legth: 180 cm,
 - Width: 65 cm made from foam,
 - Safe working load and lifting capacity of 240 kg (528 lbs),
 - Adjustable angle electric backrest from 0° to 80°,
 - Adjustable angle of footrest from 0° to -85°,
 - Excellent height range from 50 cm to 105 cm,
- **Additional Equipment (extra paid):**
 - Breathing hole and plug for prone treatments (6004),
 - Paper Roll Handle (6005),
 - Removable head/knee support (6006),
 - Removable headrest - pillow (6007),
 - Gynaecology Stirrups (6044),
 - Gynaecology leg support (6045),
 - Base cover (6052),
 - IV Pole (6054),
 - Battery Back-up (6057),
 - Double Foot switch (6207),
- Technical Passport,
- 24-month warranty for electric

