



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordynator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **Medicare Multi-Couch - Dual Footrest (SM0575) (New)**

- New,
- Made in UK,
- Electric Elevation + Electric Backrest,
- Retract Wheel Base,
- Anti-bacterial powder coated framework,
- Upper limb support,
- Vertical lifting design with stability throughout,
- TiMotion electric linear actuators,
- 90 degree dropping extendable foot sections,
- Remote Control,
- Technical Specifications:
  - Length: 180 cm,
  - Width: 65 cm made from foam,
  - Safe working load and lifting capacity of 240 kg (528 lbs),
  - Adjustable angle electric backrest from 0° to 80°,
  - Adjustable angle of footrest from 0° to -85°,
  - Excellent height range from 52 cm to 106 cm,
- **Additional Equipment (extra paid):**
  - Breathing hole and plug for prone treatments (6004),
  - Paper Roll Handle (6005),
  - Removable headrest - pillow (6007),
  - Protective Base cover (6052),
  - Battery Back-up (6057),
  - Standard operators chair (6100),
  - Premium operators chair (6150),
  - Triple foot switch (6207),
- Technical Passport,
- 24-month warranty for electric

