



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)




## Medicare Orthopaedic Couches (SM5565) (New)

- New,
- Made in UK,
- Electric,
- Retract Wheel Base,
- Anti-bacterial powder coated framework,
- Vertical lifting design with stability throughout,
- TiMotion electric linear actuators,
- 90 degree dropping extendable foot sections,
- Remote Control,
- Technical Specifications:

- Length: 180 cm or 185 cm with unfolded footrest,
- Width: 65 cm made from foam,
- Safe working load and lifting capacity of 240 kg (528 lbs),
- Adjustable angle electric backrest from 0° to 80°,
- Adjustable angle of footrest from 0° to -85°,
- Excellent height range from 52 cm to 106 cm,

### • Additional Equipment (extra paid):

- Orthopaedic arm support (5006),
- Orthopaedic full leg seat module (5007),
- Orthopaedic knee support accessory (5008),
- Orthopaedic leg support attachment (5009),
- Breathing hole and plug for prone treatments (6004),
- Removable headrest (6007),
- Paper Roll handle (6005),
- Protective Base frame shrouding (6052),
- Battery Back-up (6057),
- Operators stool (6103),
- Deluxe saddle chair (6107),
- Foot switch (6200),
- Technical Passport,
- 24-month warranty



### TAPICERKI STANDARD (W CENIE)



### TAPICERKI STANDARD + [ZA DOPŁATĄ]

