



**PRAISTON Sp. z o.o.**  
ul. Górska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## Medicare Orthopaedic Couches (SM5585) (New)


- New,
- Made in UK,
- Electric height, tilt and backrest operation,
- Basic Equipment:
  - Removable headrest (6007),
  - Paper Roll handle (6005),
- Retract Wheel Base,
- Adjustable Trendelenburg position,
- Anti-bacterial powder coated framework,
- Vertical lifting design with stability throughout,
- TiMotion electric linear actuators,
- 90 degree dropping extendable foot sections,
- Remote Control,
- Technical Specifications:
  - Length: 180 cm or 185 cm with unfolded footrest,
  - Width: 65 cm made from foam,
  - Safe working load and lifting capacity of 240 kg (528 lbs),
  - Adjustable angle electric backrest from 0° to 80°,
  - Adjustable angle of footrest from 0° to -85°,
  - Excellent height range from 52 cm to 106 cm,
- **Additional Equipment (extra paid):**
  - Orthopaedic arm support (5006),
  - Orthopaedic full leg seat module (5007),
  - Orthopaedic knee support accessory (5008),
  - Orthopaedic leg support attachment (5009),
  - Breathing hole and plug for prone treatments (6004),
  - Protective Base frame shrouding (6052),
  - Battery Back-up (6057),
  - Operators stool (6103),
  - Deluxe saddle chair (6107),
  - Triple Foot switch (6209),
- Technical Passport,
- 24-month warranty



TAPICERKI STANDARD (W CENIE)



TAPICERKI STANDARD (+ZA DOPŁATĄ)



Height range 52cm to 106cm



