



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Medicare Orthopaedic Couches (SM5585D) (New)

- New,
- Made in UK,
- Deluxe model with base shroud, orthopaedic accessories and stool,
- Basic Equipment:
 - Removable headrest (6007),
 - Paper Roll handle (6005),
 - Protective Base cover (6052),
 - Matching color operators chair (6103),
- Retract Wheel Base,
- Adjustable Trendelenburg position,
- Anti-bacterial powder coated framework,
- Vertical lifting design with stability throughout,
- TiMotion electric linear actuators,
- 90 degree dropping extendable foot sections,
- Remote Control,
- Technical Specifications:
 - Length: 180 cm or 185 cm with unfolded footrest,
 - Width: 65 cm made from foam,
 - Safe working load and lifting capacity of 240 kg (528 lbs),
 - Adjustable angle electric backrest from 0° to 80°,
 - Adjustable angle of footrest from 0° to -85°,
 - Excellent height range from 52 cm to 106 cm,
- **Additional Equipment (extra paid):**
 - Orthopaedic arm support (5006),
 - Orthopaedic full leg seat module (5007),
 - Orthopaedic knee support accessory (5008),
 - Orthopaedic leg support attachment (5009),
 - Breathing hole and plug for prone treatments (6004),
 - Battery Back-up (6057),
 - Deluxe saddle chair (6107),
 - Triple Foot switch (6209),
- Technical Passport,
- 24-month warranty



TAPICERKI STANDARD (W CENIE)



TAPICERKI STANDARD (ZA DOPŁATĄ)



