



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl


Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Phlebotomy Couch (SM9556P) (New)

- New,
- Made in UK,
- Electric with electric backrest and tilt,
- Patient handgrip supports – to aid mount/dismount from the couch and provide patient security,
- Fully adjustable phlebotomy arms for use on either side of couch,
- Fitted with paper roll holder,
- Compact frame design with retractable wheel system,
- Adjustable foot for stability on uneven surfaces,
- Purpose made precision bearings on all pivot points,
- The seat section is inclined for patient security,
- Clearance beneath frame for mobile hoist,
- TiMotion electric linear actuators,
- Technical Data:
 - Length: 180 cm,
 - Standard width 65 cm (25"),
 - Safe working load and lifting capacity (from minimum height) of 240 kg (528 lbs = 37.7 stone),
 - Adjustable angle backrest from -10 to +85 degrees positive,
 - Adjustable angle feet plates 0° to -85°,
 - Excellent height range from 50 - 102 cm (19"-39"),
- **Additional Equipment (extra paid):**
 - Removable head cushion (6007),
 - Protective base frame shrouding (6052),
 - High Operators Chair (6101),
 - Deluxe Saddle Seat (6107),
 - Standard Operators Stool (6103),
- Technical Passport,
- 24-month warranty for electrics and gas springs



Height range 52cm to 106cm



TAPICERKI STANDARD [W CENIE]



