



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl



Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Phlebotomy Couches (SM9586P) (New)

- New,
- Made in UK,
- Electric positioning for height, tilt and backrest operation,
- Adjustable angle of the footrest controlled manually by gas spring,
- Compact frame design with retractable wheel system (RWD),
- Basic equipment:
 - Fully adjustable phlebotomy arms for use on either side of couch,
 - Paper roll handle,
- Base frame powder coated with antibacterial layer,
- Adjustable for Trendelenburg position,
- Clearance beneath frame for mobile hoist,
- TiMotion electric linear actuators,
- Remote Control,
- Technical Specifications:
 - Length: 180 cm,
 - Width: 65 cm made from foam,
 - Safe working load and lifting capacity of 240 kg (528 lbs),
 - Adjustable angle electric backrest from 0° to 80°,
 - Adjustable angle footrest from 0° to - 85°,
 - Excellent height range from 52 cm to 106 cm,
- **Additional Equipment (extra paid):**
 - Short base cover (6052),
 - Battery back-up (6057),
 - High operators chair (6101),
 - Deluxe chair (6107),
 - Operators stool (6103),
 - Triple foot switch (6209),
- Technical Passport,
- 24-month warranty for electric



TAPICERKI STANDARD (W CENIE)

