




Medicare Echocardiography Couches (SM3590) (New)

- New,
- Made in UK,
- Electric positioning for height, tilt and backrest operation,
- Adjustable angle of the footrest controlled manually by gas spring,
- Adjustable for Trendelenburg position,
- Compact frame design with retractable wheel system (RWD),
- Basic equipment:
 - Removable headrest (6007),
 - Paper roll handle (6005),
- Base frame powder coated with antibacterial layer,
- TiMotion electric linear actuators,
- Remote Control,
- Technical Specifications:
 - Length: 208 cm,
 - Width: 65 cm made from foam,
 - Safe working load and lifting capacity of 240 kg (528 lbs),
 - Adjustable angle electric backrest from 0° to 80°,
 - Adjustable angle footrest from 25° to 70°,
 - Excellent height range from 50 cm to 103 cm,
- **Additional Equipment (extra paid):**
 - Battery back-up (6057),
 - Patient Lap Tray (6123),
 - Triple foot switch (6209),
- Technical Passport,
- 24-month warranty for electric



TAPICERKI STANDARD (W CENIE)

