



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Medicare Ultrasound Couches (SM9686D) (New)


- New,
- Made in UK,
- Electric positioning for height, tilt and backrest operation,
- Adjustable angle of the footrest controlled manually by gas spring,
- Adjustable for Trendelenburg position,
- Compact frame design with retractable wheel system (RWD),
- Basic equipment:
 - Removable headrest (6007),
 - Paper roll handle (6005),
 - Short base cover (6052),
 - Matching operators chair (MC6158),
- Base frame powder coated with antibacterial layer,
- TiMotion electric linear actuators,
- Remote Control,
- Technical Specifications:
 - Length: 180 cm,
 - Width: 70 cm made from foam,
 - Safe working load and lifting capacity of 240 kg (528 lbs),
 - Adjustable angle electric backrest from 0° to 80°,
 - Adjustable angle footrest from 0° to - 85°,
 - Excellent height range from 52 cm to 106 cm,
- **Additional Equipment (extra paid):**
 - Breathing hole and plus for phone treatments (6004),
 - Side support loops (6041),
 - Lithotomy stirrups (6044),
 - Leg support for gynaecology procedures (6045),
 - IV Pole (6054),
 - Battery back-up (6057),
 - Standard operators chair (6100),
 - Triple foot switch (6209),
 - Surgeons&Sonographers chair (MC6159),
- Technical Passport,
- 24-month warranty for electric



Model 6044



Model MC6158 Sonographers



TAPICERKI STANDARD (W CENIE)

